

MULTIPLE DRAWING NEEDLE

STERILE, SINGLE-USE FOR BLOOD COLLECTION



MULTIPLE DRAWING NEEDLE

STERILE, SINGLE-USE FOR BLOOD COLLECTION

- For multiple blood drawings used in combination with specific needle holders and vacuum blood collection tubes
- Three-beveled needle to minimize discomfort
- Retractable rubber sheath covering needle-end that penetrates the blood tube
- Designed to avoid kick-back at rubber sheath
- Not made with latex, DEHP or PVC
- Available in 18 G, 20 G, 21 G, 22 G
- For use by healthcare professionals only

CE 0123

Class IIa Medical Devices

Rule 6

MDD 93/42/EEC

UMDNS 12736

PRODUCT COMPLIANCE

- Complies with the following norms, directives, and regulations:
 - ISO 6009
 - ISO 7864
 - ISO 9626
 - ISO 10993-1
 - ISO 10993-4
 - ISO 10993-5
 - ISO 10993-7
 - ISO 10993-10
 - ISO 10993-11
 - ISO 11607-1
 - ISO 11607-2
 - EN ISO 13485
 - ISO 15223-1

MANUFACTURING DETAILS

Legal manufacturer: Nipro Corporation

Country of origin: JAPAN

STERILIZATION AND SHELF LIFE

Sterilization

ETO (Ethylene oxide)



Single-use only



Indicator sticker


Each outer box (or shipping carton) contains a chemical indicator that indicates sterility. This indicator comes in the form of a blue sticker that turns red when sterilized.

Shelf life

3 years



MATERIALS USED



PRODUCT RANGE OVERVIEW

Reference code	Description	Color
NM- 1825-EJ	Multi Draw Ndl 18Gx1" (1,2x25mm)	Pink
NM-1838-EJ	Multi Draw Ndl 18Gx1 1/2" (1,2x38mm)	Pink
NM -2025-EJ	Multi Draw Ndl 20Gx1" (0,9x25mm)	Yellow
NM-2038-EJ	Multi Draw Ndl 20Gx1 1/2" (0,9x38mm)	Yellow
NM-2125-EJ	Multi Draw Ndl 21Gx1" (0,8x25mm)	Green
NM-2138-EJ	Multi Draw Ndl 21Gx1 1/2" (0,8x38mm)	Green
NM-2225-EJ	Multi Draw Ndl 22Gx1" (0,7x25mm)	Black
NM-2238-EJ	Multi Draw Ndl 22Gx1 1/2" (0,7x38mm)	Black

PACKAGING DETAILS

Reference code	Inner box				Intermediate box			Outer box (mm)					
	Quantity (pcs)	L (cm)	W (cm)	H (cm)	L (cm)	W (cm)	H (cm)	Quantity (pcs)	L (cm)	W (cm)	H (cm)	Weight (kg)	M ³
NM-1825-EJ	100	8,1	16,4	5,2	42,3	22,4	17,1	6000	54	44,5	24,5	16	0,058
NM-1838-EJ	100	8,1	16,4	5,2	42,3	22,4	17,1	6000	54	44,5	24,5	16	0,058
NM-2025-EJ	100	8,1	16,4	5,2	42,3	22,4	17,1	6000	54	44,5	24,5	15,54	0,058
NM-2038-EJ	100	8,1	16,4	5,2	42,3	22,4	17,1	6000	54	44,5	24,5	15,7	0,058
NM-2125-EJ	100	8,1	16,4	5,2	42,3	22,4	17,1	6000	54	44,5	24,5	15,32	0,058
NM-2138-EJ	100	8,1	16,4	5,2	42,3	22,4	17,1	6000	54	44,5	24,5	15,50	0,058
NM-2225-EJ	100	8,1	16,4	5,2	42,3	22,4	17,1	6000	54	44,5	24,5	15,18	0,058
NM-2238-EJ	100	8,1	16,4	5,2	42,3	22,4	17,1	6000	54	44,5	24,5	15,61	0,058

Languages

12 languages

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Finnish (FI), Swedish (SV), Danish (DA), Norwegian (NO), Greek (EL).










Transport conditions

Closed and dry

Storage conditions

Open the packaging only immediately before use to guarantee sterility.

Label details

	Individual unit	Inner box	Intermediate box	Outer box
Product name	✓	✓	✓	✓
Product size	✓	✓	✓	✓
Quantity (pieces)	—	✓	✓	✓
	✓	✓	✓	✓
	✓	—	—	—
	✓	✓	✓	✓
	✓	✓	✓	✓
	—	✓	✓	✓
	—	✓	✓	✓
	—	✓	✓	✓
	—	✓	✓	✓
	—	✓	✓	✓
	—	✓	✓	✓

